

TECH TIPS

CARDIAC DEVICES

UPDATED LABELING

FOR 3T

Philips Specific

B_{1+RMS} as a Condition of Use

Because B_{1+RMS} is now recommended as a supplemental metric to SAR, Medtronic has begun to have device labeling approved with B_{1+RMS} limits.

For a full list of devices and leads approved for the MRI environment, refer to our MR-conditional Cardiac Device Summary Chart.

SureScan™ cardiac device 3T labeling is a little different than 1.5T labeling in that it **uses B_{1+RMS} rather than SAR when landmarking/centering below C7**. When scanning these devices at 1.5 T, SAR must be limited to the Normal Operating Mode. It is important to remember that the B_{1+RMS} limit applies to 3T ONLY. When utilizing B_{1+RMS} as a condition of use at 3T, one can utilize either Normal or First Level Controlled Mode for SAR. However, the displayed B_{1+RMS} value for each scan is not to exceed 2.8 μT (micro Tesla) when landmarking/centering below C7. When landmarking/centering above C7, either Normal or First Level Controlled mode may be used for SAR and there are no restrictions regarding B_{1+RMS} . If your 3T MR system does not display B_{1+RMS} , then only studies may be performed which are landmarked/centered above C7. **Centering/Landmarking and/or scanning below C7 is not allowed for 3T systems that do not display B_{1+RMS} .**

At 3T, most of your current protocols will need little to no modification to comply with the 2.8 μT B_{1+RMS} condition of use. As a rule of thumb, the higher the estimated SAR, the higher the likelihood that some modifications may need to be made. However at 3T, the likelihood that a sequence would need modification to comply with the 2.8 μT limit is still very low.

It should also be noted that there are no restrictions on the use of local transmit/receive coils for imaging of the head or the extremities, which includes no B_{1+RMS} restrictions. For example, if one is using a knee transmit/receive coil with a sequence requiring > 2.8 μT , this would be acceptable since the whole body coil is not being used as the transmit coil.



How one modifies scan parameters to affect the B_{1+RMS} may vary slightly between MR system brands based on available options. However, in general, whatever you do to reduce SAR on your system will likely reduce the B_{1+RMS} . If your system displays B_{1+RMS} , you will find it in the same area in which you find SAR.

Philips 3T MR Systems

On the Philips system, B_{1+RMS} is displayed alongside the estimated SAR values as shown in Figure 1. In particular it is the "Max B_{1+RMS} " that is to be monitored.




Figure 1

As previously mentioned, most parameters that affect SAR will affect the B_{1+RMS} . Therefore, setting the **SAR mode to "low"** will reduce the B_{1+RMS} . However, it is highly unlikely this will be necessary to stay below 2.8 μT . Once you have adjusted the parameters to obtain a B_{1+RMS} value of 2.8 μT or less, you can save that sequence in your protocols. Unlike SAR, the B_{1+RMS} value will be the same the next time you recall that sequence (provided none of the parameters are altered from when it was saved, including the number of slices).

Specific Examples for Modifying B_{1+RMS} on Philips 3T Systems

Reducing the number of slices without increasing the TR or increasing the TR without increasing the number of slices will reduce the B_{1+RMS} . Reducing the number of slices is often not practical so increasing the TR is the more likely choice between the two. Depending on other parameters, you may have to significantly increase the TR. Besides the impact on scan time, significant increases in TR are also not practical when T1-Weighted spin echo sequences are desired.

Fast spin echo series consist of a 90-degree pulse followed by a "train" of refocusing pulses, generally said to be 180-degree pulses. In reality these are rarely 180-degree pulses but rather 170-degrees or even less. The number of echoes generated by the refocusing pulses is referred to as the **TSE factor**. To reduce the B_{1+RMS} you can reduce the TSE factor (leaving all other parameters unchanged).

| Summary | Geometry | Contrast | Motion | Dyn/Ang | Postproc | Offc/Ang | Coils | Conflicts | << |
|--------------------|----------|--------------|--------|---------|---------------------------|----------|---------------------|-----------|----|
| ultrashort | | no | | | Act. TE (ms) | | 80 | | |
| fid reduction | | default | | | ACQ matrix M x P | | 360 x 256 | | |
| Echoes | | 1 | | | ACQ voxel MPS (mm) | | 0.50 / 0.70 / 5.00 | | |
| partial echo | | no | | | REC voxel MPS (mm) | | 0.47 / 0.47 / 5.00 | | |
| TE | | user defined | | | Scan percentage (%) | | 71.3 | | |
| (ms) | | 80 | | | Packages | | 6 | | |
| Flip angle (deg) | | 90 | | | Min. slice gap (mm) | | 5 | | |
| Refocusing control | | yes | | | Optimal slices | | 12 | | |
| angle (deg) | | 120 | | | WFS (pix) / BW (Hz) | | 1.370 / 317.1 | | |
| TR | | range | | | Full flow comp. | | yes | | |
| minimum (ms) | | 3000 | | | TSE es / shot (ms) | | 8.9 / 151 | | |
| maximum (ms) | | 5000 | | | Min. TR (ms) | | 3966 | | |
| Halfscan | | no | | | SAR / local torso | | < 50 % | | |
| Water-fat shift | | user defined | | | Whole body / level | | < 1.3 W/kg / normal | | |
| (pixels) | | 1.8 | | | SED | | < 0.9 kJ/kg | | |
| RF Shims | | adaptive | | | B1+rms / Coil Power | | 1.42 uT / 37 % | | |
| Shim | | auto | | | Max B1+rms | | 1.43 uT | | |
| mDIXON | | no | | | PNS / level | | 59 % / normal | | |
| Fat suppression | | no | | | dB/dt | | 35.2 T/s | | |
| Water suppression | | no | | | Sound Pressure Level (dB) | | 23.2 | | |

Figure 2

You also have the option to alter the flip angle of the refocusing pulse. (Figure 2). You will find this option under the contrast tab. Selecting "yes" for Refocusing control will open up a selection to allow you to reduce the angle of the refocusing pulse. Care should be taken when reducing the refocusing angle. If it is reduced significantly, you may see a reduction in SNR and/or and alteration of image contrast. Experimentation may be the best way to determine the lower limit for any given sequence. As you reduce either the TSE factor or the refocusing flip angle, the B_{1+RMS} will be reduced.

Please note that as of 2016, all Philips 3T systems specify a maximum B_{1+RMS} of 2.3 μT in their compatibility technical data sheet, thus any modifications of any kind are not expected to be needed.

Summary

- 3T labeling uses B_{1+RMS} as the condition of use relative to RF power when landmarking/centering below C7.
- When **landmarking/centering above C7**, there are no restrictions for B_{1+RMS} and either Normal or First Level Controlled mode for SAR may be selected.
- When **landmarking/centering below C7**, the displayed B_{1+RMS} value should be less than or equal to 2.8 μT (micro-Tesla). Either Normal or First Level Controlled SAR mode may be selected.
- In the event your 3T system software does not display B_{1+RMS} , only studies in which the landmark is above C7 may be performed on a 3T system.
- The use of B_{1+RMS} as a metric for RF heating provides for greater flexibility in pulse sequence and parameter selection.
- Most parameters which affect SAR will affect B_{1+RMS} .
- Once a sequence has been modified to have a B_{1+RMS} value of 2.8 μT or less, it can be saved in the site's protocols.
- As long as the parameters affecting B_{1+RMS} are not modified, sequences saved with a specific B_{1+RMS} value will remain unchanged patient-to-patient.

Brief Statement

Medtronic SureScan™ Portfolio for 1.5T and 3T MR-conditional Use

Medtronic SureScan products and systems are MR Conditional, and as such are designed to allow patients to undergo MRI under the specified conditions for use.

Pacing, ICD, and CRT-D Systems: When programmed to On, the MRI SureScan feature allows the patient to be safely scanned while the device continues to provide appropriate pacing. A complete transvenous SureScan system, which is a SureScan device with appropriate SureScan lead(s), is required for use in the MR environment. For ICD and CRT-D Systems, when a single coil SureScan defibrillation lead is used, a Medtronic DF-1 pin plug must be secured in the SVC port to make a complete SureScan DF-1 defibrillation system. To verify that components are part of a SureScan system, visit <http://www.mrisurescan.com/>. Any other combination may result in a hazard to the patient during an MRI scan.

See the MRI SureScan Technical Manual before performing an MRI Scan and Device Manual for detailed information regarding the implant procedure, indications, contraindications, warnings, precautions, and potential complications/adverse events. For further information consult Medtronic's website at www.medtronic.com or www.mrisurescan.com.

Products that appear on this website may not all be approved in your country. Please contact your local affiliate for further information.

Medtronic

99 Hereford Street
Brampton, ON L6Y 0R3
Canada

Tel: (905) 460-3800
Fax: (905) 826-6620

medtronic.ca

UC201700631b EC ©2017
Medtronic. Minneapolis, MN. All Rights
Reserved. Printed in USA. 06/2017

Medtronic