

MR-CONDITIONAL CARDIAC DEVICE SUMMARY CHART

January 2018

- The chart below contains all Medtronic cardiac devices Health Canada licensed for MRI scans under specific conditions for use
- If a model number and/or lead length is not listed, then it is not Health Canada licensed for the MR environment

All devices listed are licensed for 1.5T and 3.0T scanning.

THERAPY	PRODUCT	MODEL NUMBER	MR CONDITIONAL 	SURESCAN LEADS
PACEMAKERS	Advisa MRI™ SureScan	A3DR01 A3SR01	Yes—If complete system is implanted with a SureScan™ pacemaker and SureScan lead(s)	Pacing leads: 5076: 35, 45, 52, 58, 65, 85 cm 4076: 35, 45, 52, 58, 65, 85 cm 5086MRI: 45, 52, 58 cm 3830: 59, 69, 74 cm 4074: 52, 58cm 4574: 45,53 cm
	Ensura MRI™ SureScan	EN1DR01 EN1SR01		
	Micra TPS	MC1VR01		
CARDIAC DEFIBRILLATORS	Evera MRI™ XT	DDMB2D4 DMB2D1 DVMB2D4 DVMB2D1	Yes—If complete system is implanted with a SureScan ICD and SureScan lead(s)	Defibrillation leads: 6935, 6947, 6935M, and 6947M Lengths 58/65cm for 6935 and 6947, and 55/62cm for 6935M and 6947M
	Evera MRI S	DDMC3D4 DDMC3D1 DVMC3D4 DVMC3D1		
	Visia AF MRI XT	DVFB2D4		
	Visia AF MRI S	DVFC3D4		
INSERTABLE CARDIAC MONITORS	Reveal™ XT	9529	Yes	Not Applicable
	Reveal™ DX	9528	Yes	
	Reveal LINQ™	LNQ11	Yes	Not Applicable
CARDIAC RESYNCHRONIZATION THERAPY DEFIBRILLATORS	Claria MRI™ Quad CRT-D	DTMA2QQ	Yes—If complete system is implanted with a SureScan CRT-D and SureScan lead(s) or the Model 6725 pin plug	CRT leads: 4196, 4296 and 4396 Lengths: 78 and 88 cm 4298, 4398 and 4598 Lengths: 78 and 88 cm See above for defibrillation and pacing leads.
	Claria MRI™ CRT-D	DTMA2D4		
	Amplia MRI™ Quad CRT-D Amplia MRI™ CRT-D	DTMB2QQ DTMB2D4		
	Compia MRI Quad CRT-D Compia MRI CRT-D	DTMC2QQ DTMC2D4		

To obtain the MR Conditions for Use, go to www.manuals.medtronic.com or www.mrisurescan.com/productsearch.

Brief Statement

SureScan™ Pacing Systems and Defibrillation

Medtronic SureScan pacing systems are MR-Conditional and as such are designed to allow patients to undergo MRI under the specified conditions for use. A complete SureScan pacing system, which consists of an approved combination (see <http://www.mrisurescan.com/ca>) MRI SureScan device with SureScan lead(s), is required for use in the MRI environment. Consult the device manuals to ensure all system components are MR-Conditional. (Some products may not be licensed in accordance with Canadian law.)

See the MRI SureScan Technical Manual before performing an MRI Scan and Device Manual for detailed information regarding the implant procedure, indications, contraindications, warnings, precautions, and potential complications/adverse events.

Medtronic of Canada
99 Hereford Street
Brampton, ON L6Y 0R3
Canada
Tel: (905) 460-3800
Fax: (905) 826-6620

Toll-free: 1 (800) 505-4636
(24-hour technical support for
physicians and medical professionals)

UC201204952ke EC©Medtronic
2016. Minneapolis, MN. All Rights
Reserved.

medtronic.ca

Medtronic